







DIANA^{ADV}

RE-CREATE



DIANA^{ADV}

SUPER PULSE
Co2 FRACTIONAL



Having laser output mode as super pulse, **Diana - ADV** has its applications in skin renewal and pigmentation removal. The product is an effective method to get rid of the age spots, blemishes, speckles and warts. Its scope of treatment also lies in treating skin scars and vaginal problems. The product is well-equipped with gynaecology laser treatment to sort out the issues of vaginal tightening, urine incontinence and vaginal whitening.

TECHNICAL SPECIFICATIONS


LASER WAVELENGTH	: Regular mode 30/40/60W. Fractional mode 0.1-300mj
OUTPUT POWER	: 10600nm
LASER SOURCE	: RF Excited Tube
LASER OUTPUT MODE	: Super pulse, CW Single pulse, Repeated pulse
WORKING MODE	: Fractional mode, Surgical mode, Vaginal mode (optional)
SCANNING DIRECTION	: Sequence, Divergence, Random scanning
SCANNING SIZE	: Up to 20*20mm
SPOT DIAMETER	: 0.08-0.12mm
PULSE WIDTH	: Fractional mode: 0.1-10ms, Surgical mode: 05-1000ms
AIMING BEAM	: 650nm red diode laser, laser power 0.1-5mW
OPTICAL TRANSMISSION SYSTEM	: 7 joints, 360o rotation
COOLING SYSTEM	: Air cooling
INPUT POWER	: < 600W



GYNAECOLOGY LASER TREATMENT


KEY FEATURES :

✓ 7 VARIABLE TREATMENT GRAPHICS, ADJUSTABLE SHAPES, SIZES AND SPACING



TREATMENT HEADS FOR 3 SYSTEMS: NORMAL, FRACTIONAL AND GYNECOLOGY

BEFORE & AFTER :



DIAN ADV

 #GF, KH no 71/4, Vill. Bamnoli, Dwarka Sec 28, Near By Safe Shop
New Delhi-110071

 Phone : 011-28041611

 info@woxmen.com, woxmenlaser@woxmen.com

 www.woxmen.com

Regional Office :

North India - Mob.: +91-9999416987, 8800835930

South India - Mob.: +91-7892485400

East India - Mob.: +91-8209520887, 9166608900

West India - Mob.: +91-9611101560

