








YOUR BEAUTY, OUR DUTY.





IRIS FV

PICO LASER

IRIS - FV is a one-stop solution for the removal of all kinds of marks on your body. The device consists of 7 joined arms with a motor-driven auto-switching facility to treat dermis pigmentation, man-made pigmentation, vascular lentiginos and other marks. Its laser technology has its application in skin rejuvenation. Being user-friendly and easily manageable, the product is guaranteed to perform efficiently without any side effects. The best fact of the product is the tiny sense of pain due to which Anaesthesia is not required in most cases.

TECHNICAL SPECIFICATIONS

WAVELENGTH	: 1064nm , 532nm
OUTPUT ENERGY	: 1064nm: 100-1500mj, 532nm: 100-800mj
PULSE DURATION/WIDTH	: 8(+/-) 12ns / 6ns / 600ps
REPETITION RATE	: 1-10 Hz Adjustable
SPOT SIZE	: 1-10 mm Adjustable
LIGHT GUIDE BEAM	: 655nm red laser
COOLING SYSTEM	: Close Cycle Water+ Air Heat Exchanger
LASER DELIVERY	: 7 Jointed Articulated Arm
CONTROLLING SYSTEM	: 10.4 inch TFT color touch screen
POWER REQUIREMENT	: AC:220V/110V±10%, 50Hz/60Hz
LANGUAGE OPTIONS	: English / Korean
MACHINE DIMENSIONS	: 680mm X 400mmX820mm
NET WEIGHT	: 75Kgs
WARRANTY	: 2 Year from Installation



IRIS FV



KEY FEATURES :

- ✓ 7 joint articulated arm helps treating patient easily and reduces overall treatment time
- ✓ Red diode laser beam helps targeting desired area with pin point accuracy
- ✓ Motor driven auto switching for 1064nm and 532nm wavelength
- ✓ Adjustable spot size through focus adjustment lens
- ✓ High power cooling system for continuous operations
- ✓ No harm to hair follicle and normal skin, leaves no scar after treatment
- ✓ High quality controlled solid state laser in compliance with international standards
- ✓ No downtime and interruption of routine activities, no side effects
- ✓ User friendly interface and easy to operate
- ✓ Tiny sense of pain, No Anesthesia required in most of cases




7 Joint
Articulated Arm



Japanese
Cooling System



Red Diode
Target Beam



Spot Size
Adjustment

SCOPE OF TREATMENT

- ➡ Dermis Pigmentation – Birth mark such as : Ota nevus, Ito nevus, black nevus, blue nevus
- ➡ Treatment of sun-burn spot, age spot, speckle, freckle, senior fleck
- ➡ Treatment of Manmade Pigment: Blue, Black, Red, Brown Tatoo, Lip line, Eyebrow line
- ➡ Epidermis Pigmentation, Speckles: Freckles, Senior Freckles, Coffee Spot
- ➡ Vascular Lentiginous: Vascular lesion, Spider Vessel, Spider Veins
- ➡ Skin Rejuvenation with laser peeling

 #GF, KH no 71/4, Vill. Bamnoli, Dwarka Sec 28, Near By Safe Shop
New Delhi-110071

 Phone : 011-28041611

 info@woxmen.com, woxmenlaser@woxmen.com

 www.woxmen.com

Regional Office :

North India - Mob.: +91-9999416987, 8800835930

South India - Mob.: +91-7892485400

East India - Mob.: +91-8209520887, 9166608900

West India - Mob.: +91-9611101560

