



Pelviwave

EMS Chair For Pelvic Muscles



PelviWave

EMS Chair For Pelvic Muscles

The PelviWave chair is a non-invasive and painless device for the treatment of pelvic floor muscle weakness, urinary incontinence and vaginal muscle laxity, utilizing High Intensity Focused Electro- Magnetic Technology to stimulate deep pelvic floor muscles and restore neuromuscular stimulation!

Key effectiveness is based on focused electromagnetic energy, in-depth penetration and stimulation of the entire pelvic floor area.

A single PelviWave chair session brings more than 10,000 supramaximal pelvic floor muscle contractions compared to common Kegel exercises, which are extremely important in muscle re-education of incontinence patients.


How Does It Work?

For Women




Before

Pelvic Floor Muscles Insufficient Support Pelvic Organs And Affect Bladder Control.



Treatment


Using High Intensity Focused Electro-Magnetic Technology Effectively Stimulates Pelvic Floor Muscles with Thousands of Supramaximal Contractions.



After


Stimulation Leads To Regained Control Over Pelvic Floor Muscles And Bladder.

For Men



Before

Pelvic Floor Muscles are Insufficiently Support Pelvic Organs and Affect Bladder Control.



Treatment

Pelviwave Effectively Stimulates The Pelvic Floor Muscles Using High Intensity Focused Electro-Magnetic Technology



After

Stimulation Leads To Regain Control Over Pelvic Floor Muscles and Bladder

Advantages



Non-invasive and no recovery time



FDA Cleared Treatment For Incontinence



75% Pad Reduction



Fully clothed throughout treatment



Walk In And Walk Out



Comfortable & Painless

Technical Specifications

Power	3000W
Output Method	Pulse Output
Control Method	Touch Control
Screen Size	7-inch Capacitive Screen
Output Frequency	1-150HZ
Net Weight	32KG
Dimension	560x620x860MM
Package Size	710x630x670MM
Package Weight	45KG
Voltage	220V/110V

LEATHER BACKREST
High-end and comfortable


7-INCH REMOVABLE TOUCH SCREEN
The screen can be held directly in your hand, which is very convenient to use

STRONG MEDICAL WHEELS
The maximum load-bearing capacity can reach 200kg


SUPER LARGE MAGNETIC COIL
30cm large magnetic coil, maximum output magnetic force can reach 30uH

3000W STRONG POWER SUPPLY
The unique 3000W ultra-high power supply ensures high-intensity output of the machine and the best treatment effect

Treatment Results




Before



Before



Atrophic, Postpartum
Pelvic Floor muscle



Normal, Healthy
Pelvic Floor Muscle